

US-PAT-NO: 6093530

DOCUMENT-IDENTIFIER: US 6093530 A

TITLE: Non-calcific biomaterial by glutaraldehyde followed
by oxidative
fixation

----- KWIC -----

The present invention relates to a method of preparing biomaterial for prosthetic use that has excellent mechanical properties and high calcification resistance. The method includes the processing of tissues with glutaraldehyde, by standard methods, to fix the tissue. The method then calls for oxidizing the tissue using photosensitive dyes or other means. The resulting tissue has the favorable mechanical properties associated with glutaraldehyde fixation and the resistance to calcification associated with oxidative fixation.

6. The process of claim 1 wherein the fixation process comprises exposing the biomaterial to glutaraldehyde in phosphate buffered saline solution for about 12 to 36 hours at a temperature from about 1.degree. to 10.degree. C., followed by exposing the biomaterial to glutaraldehyde in phosphate buffered saline solution for about 12 to 36 hours at a temperature from about 15.degree. to 30.degree. C., followed by exposing the biomaterial to glutaraldehyde in phosphate buffered saline solution for about 12 to 36 hours at a temperature from about 30.degree. to 44.degree. C., and the photooxidation process comprises exposing the biomaterial to a solution comprising dissolved oxygen and one or more of a compound selected from methylene blue,

methylene green,
rose bengal, riboflavin, proflavin, fluorescein, rosin, and
pyridoxal-5-phosphate, and then exposing the tissue in the
solution to light in
quantities sufficient to cause crosslinking of the tissue.

US-PAT-NO: 5746775

DOCUMENT-IDENTIFIER: US 5746775 A

TITLE: Method of making calcification-resistant
bioprosthetic tissue

----- KWIC -----

The overall protein composition and valvular morphology of porcine aortic valves are unaffected by alcoholic treatment as demonstrated by complete amino acid analysis and electron spectroscopy for chemical analyses (ESCA). In fact, alcohol treatment enhances surface smoothing and anisotropy of porcine aortic valve leaflets resulting in a surface chemistry which is comparable to fresh leaflets. In contrast, glutaraldehyde-pretreated (control) or detergent (SDS) treated tissue show significant differences. Table 5 hereinbelow presents ESCA data of the surface carbon (C1s), nitrogen (N1s), and oxygen (O1s) concentrations (%) in porcine aortic valve specimens immersed for 24 hours in the indicated solution.

US-PAT-NO: 5934283

DOCUMENT-IDENTIFIER: US 5934283 A

TITLE: Pubovaginal sling device

----- KWIC -----

The cross-linking of the tissues used in the present invention is preferably accomplished using one or more of the following treatment agents:
glutaraldehyde, dialdehyde, glutaraldehyde starch, dialdehyde starch, an epoxy compound or ionizing radiation. Certain processes (such as heat, radiation or pH change) or agents such as halogens, enzymes, organic solvents, detergents, sodium hydroxide, hydrochloric acid, sodium hypochlorite or hydrogen peroxide) may be used to inactivate viruses with and without protein coats or to destroy BSE agent infectivity during the manufacturing process. The tissue may also be treated with a highly volatile chemical such as for example propylene oxide, to assist with the sterilization of the tissue. Sterilization may be accomplished using one or more of the following treatments:
glutaraldehyde, alcohol, propylene oxide or irradiation sterilization. The treatment of the tissue, with a combination of these materials and processes, can both cross-link the tissue and render the tissue sterile for implantation.

US-PAT-NO: 5819748

DOCUMENT-IDENTIFIER: US 5819748 A

TITLE: Implant for use in bone surgery

----- KWIC -----

There are significant differences between collagen fibres at different conditions of maturity. Where the connective tissue is in an active phase of fibrillogenesis, for example during growth or wound healing, collagen fractions can be isolated with different properties. The first fraction is extractable by neutral solutions (neutral-soluble collagen); this consists of recently synthesised tropocollagen molecules which are not aggregated or are only beginning to aggregate. The second fraction is extractable by a sodium citrate solution at pH 3.0, and is thus termed the acid-soluble collagen fraction. The third fraction found in older tissues is the insoluble fraction and can only be extracted by very vigorous methods. One basis for the difference between these fractions lies in the degree of cross-linking by oxidation to produce peroxide bridges. Collagen can also be cross-linked chemically via free amino groups, using aldehydes such as formaldehyde or glutaraldehyde or isocyanates such as hexamethylene diisocyanate. By such cross-linking, animal collagen fractions lose their antigenicity almost completely. Cross-linking of collagen fibrils in this way is for example, made use of in heart replacement surgery, where animal, e.g. porcine, valves are conditioned with glutaraldehyde for use as human pulmonary or mitral valve replacements.

US-PAT-NO: 5697972

DOCUMENT-IDENTIFIER: US 5697972 A

TITLE: Bioprosthetic heart valves having high
calcification resistance

----- KWIC -----

Carpentier teaches that inflammatory reactions including immune reactions and the degeneration of collagen and elastin are the major factors of malfunction in bioprosthetic heart valves after transplantation, and suggests a new method to eliminate these factors for clinical use [A. Carpentier, Biological Tissue in Heart Valve Replacement, M. I. Ionescu et al. (Eds). Butterworth, London, 1972]. He developed a method for treating the bioprosthetic heart valve to inhibit any inflammatory reactions by host cells, to maintain or enhance its strength and flexibility under sterilized conditions, to prevent a degeneration of collagen and elastin, and to prevent the intrusion of host cells in the transplanted valve. This method comprises taking a porcine valve aseptically, washing it with Hanks solution to remove soluble antigenic material, and then oxidizing mucopolysaccharide and glycoprotein with sodium metaperiodate to form an aldehyde group on the side chains thereof. The aldehyde group thus formed may be bound with the adjacent amine group to form an intermolecular cross-linkage. The remaining sodium metaperiodate is neutralized with ethylene glycol and then residual amine groups of the glycoprotein molecule are cross-linked with glutaraldehyde buffer solution. Finally, these cross-linkages are stabilized by a reduction with sodium

borohydride. However, this method has disadvantages that the durability of the transplanted valve is decreased by calcification during the long-term use.

A bioprosthetic tissue of a bovine pericardium was immersed in 200 ml of Hanks solution (10%) containing H.sub.2 N--PEO400--SO.sub.3 at 4.degree. C. for 2 hours, and reacted with EDC without washing. Then, the tissue in 200 ml of a 1% H.sub.2 N--PEO400--SO.sub.3 NaIO.sub.4 solution was oxidized at 4.degree. C. for 24 hours, and neutralized at 4.degree. C. for 1 hour in 200 ml of a 10% H.sub.2 N--PEO400--SO.sub.3 /ethylene glycol solution. The tissue thus treated was fixed with a glutaraldehyde solution at 4.degree. C. for a week and reduced with NaBH.sub.4 to prepare BT--PEO400--SO.sub.3.

A bioprosthetic tissue of a bovine pericardium was immersed in Hanks solution for 2 hours, and then in 50 ml of 2% glutaraldehyde at 4.degree. C. for 24 hours, and washed with a PBS solution. The tissue was added to 50 ml of a 5% H.sub.2 N--PEO2000--SO.sub.3 solution having pH 11, and the mixture allowed to react at room temperature for 24 hours. The tissue was washed with a PBS solution, and added to 200 ml of a 1% sodium metaperiodate (NaIO.sub.4) solution containing a 1:2 mixture of a 3% NaIO.sub.4 solution and Hanks solution. The tissue was added to 200 ml of a 1% NaIO.sub.4 solution, and oxidized at 4.degree. C. for 24 hours in the dark. The tissue was then washed with a PBS solution, and added to 200 ml of a 1% ethylene glycol solution to neutralize at 4.degree. C. for 1 hour. The tissue was washed with a PBS solution, and treated with an acetic anhydride solution to block the unreacted amine group. The tissue was further reacted with H.sub.2 N--PEO2000--SO.sub.3 using EDC. The tissue thus treated was fixed with a

glutaraldehyde solution
for a week, and reduced with NaBH.sub.4 to prepare
BT--PEO2000--SO.sub.3.

A bioprosthetic tissue of a porcine aortic valve was
immersed in Hanks solution
for 2 hours, and then in 50 ml of 2% glutaraldehyde at
4.degree. C. for 24
hours, and washed with a PBS solution. The tissue was
added to 50 ml of a 5%
H.sub.2 N--PEO1000--SO.sub.3 solution having pH 11, and
allowed to react at
room temperature for 24 hours. The tissue was washed with
a PBS solution, and
treated with an acetic anhydride solution to block the
unreacted amine group.
Then, the tissue was further reacted with H.sub.2
N--PEO1000--SO.sub.3 using
EDC. The tissue was washed with a PBS solution, added to
200 ml of a 1%
NaIO.sub.4 solution containing a 1:2 mixture of a 3%
NaIO.sub.4 solution and
Hanks solution. The tissue was oxidized at 4.degree. C.
for 24 hours in the
dark. The tissue thus treated was washed with a PBS
solution, and added to 200
ml of a 1% ethylene glycol solution to neutralize at
4.degree. C. for 1 hours.
Then, the tissue was fixed with glutaraldehyde solution for
a week, and reduced
with NaBH.sub.4 to prepare BT--PEO1000--SO.sub.3.